




LOMBOSTAT ΑΣΦΑΛΕΙΑΣ

Ζώνη οσφύος ελαστική, αυτοκόλλητη με επιπλέον αγκράφες μπροστά.

Συγκρατεί & στηρίζει τη μέση με 6 μπανέλες και ιμάντες χιαστί.

Οσφυαλγία, ισχιαλγία, σπονδυλοαρθροπάθεια.

ΜΕΓΕΘΟΛΟΓΙΟ



◀ Μετρήστε γύρω από τη μέση
στο ύψος του ομφαλού.

Μεγέθη	80	85	90	95	100	105	110	115
Μέση	80	85	90	95	100	105	110	115

► Για οποιαδήποτε περαιτέρω πληροφορία ή διευκρίνιση,
παρακαλώ μη διστάσετε να επικοινωνήσετε μαζί μας
στο **211 10 24 100**



Ref.no. **11300**

Χρώμα: Μπεζ

Ύψος: 22 cm

ΟΡΘΟ ΔΥΝΑΜΙΚΗ
αβε

Ναυπλίου 45 - 144 52
Μεταμόρφωση - ΑΤΤΙΚΗΣ

211 10 24 100



info@orthodynamic.gr
www.orthodynamic.gr






LOMBOSTAT with Stays

Hard tension elastic material with six metallic splints and crossed straps in the back. With velcro closure and extra claps in front for better safety.

Indications chronic lumbago, lumbago sciatica.



SIZE GUIDE



◀ Measure circumference around the waist immediately below navel

Sizes	80	85	90	95	100	105	110	115
Waist	80	85	90	95	100	105	110	115

► If you have any questions or would you like more information, please contact with us **+30 211 10 24 100**

Ref.no. **1 1300**

Color: Beige

Height: 22 cm



ORTHO DYNAMIC
sa

45 Nafpliou Str.,GR -144 52
Metamorfofi - Greece

+30 211 10 24 100



info@orthodynamic.gr
www.orthodynamic.gr

